



Estimating the actual dose delivered by intravascular coronary brachytherapy using geometrically correct 3-D modeling

Andreas Wahle, John J. Lopez,* Edward C. Pennington, Sanford L. Meeks, Kathleen C. Braddy, Theresa M. H. Brennan, John M. Buatti, James D. Rossen, Milan Sonka


The University of Iowa, Departments of Electrical and Computer Engineering, Internal Medicine, and Radiation Oncology

*The University of Chicago, Department of Medicine


Motivation

- Coronary arteries tend to develop new plaque accumulation after balloon angioplasty and stenting (*in-stent restenosis*)
- Irradiating the restenosed vessel segment decrease recurrent restenosis rate by 40-60%
- Common dose models for intravascular brachytherapy do not consider vessel curvature and eccentricity of catheter


In-Stent Restenosis:




(a)




(b)



(c)



(d)



(e)


In-Stent Restenosis (angio)





RAO





LAO



In-Stent Restenosis (IVUS)




Angio- graphy and IVUS



Position of the IVUS Catheter within Vessel


**Consequences for
Brachytherapy?**




Questions:

- What is the impact of the vessel *curvature* on dose distribution?
- What is the impact of catheter *eccentricity* on dose distribution?


Dose Distribution in 0.5-mm Layers




Simulation: 60-180° Torus




Dose Distribution 60° Torus




Dose Distribution 120° Torus



Dose Distribution 180° Torus




Results 120° Torus



Prescribed dose: 18.4 Gy

Results 120° Torus



Eccentricity of the Catheter




Scenario Inner Curvature

Scenario Outer Curvature

40mm 16-Source Train within 60mm Simulated Vessel

120-degree Curved Vessel with Eccentric Catheter




Simulation Results:

- Vessel *curvature* increases the doses delivered to the inner bend
- Catheter *eccentricity* biases dose distribution towards closer wall
- Effects may partially offset each other if catheter is at outer bend

Questions (2):

- Impact *in-vivo*?
- How to generate an accurate 3-D model of the vessel segment?
- Angiography-IVUS Fusion

Principle of Data Fusion





IVUS

Angiography

Fusion Outline


- Matching of IVUS frames on 3-D path
 - Constant pullback speed
- Determination of absolute orientation
 - Differential geometry
 - Using IVUS catheter as landmark
- Mapping of pixel and contour data

In-Vivo Example







35LAO




Finite-Element Mesh



Segment 2



Segment 1




Simulation




2 Segments
30mm train

Doses


1. Adventitia
2. Lumen



Visualization




**Geometrically
Correct Model**





Segment 2

Segment 1

**Simplified
Tubular
Model**




 Lumen
 Plaque+
Media

Geometrically Correct Model



Simplified Tubular Model



In-Vivo Results

Layer	GCM [Gy]	STM [Gy]	Difference means [%GCM]			
	(mean±SD)	(mean±SD)	All	LAD	LCX	RCA
	<i>n</i> = 10	<i>n</i> = 10	<i>n</i> = 10	<i>n</i> = 3	<i>n</i> = 2	<i>n</i> = 5
Lumen/Plaque	34.70 ± 11.93	31.66 ± 5.11	8.76	6.54	7.62	10.44
Intermediate 1	26.71 ± 8.76	24.85 ± 4.14	6.96	4.77	6.58	8.38
Intermediate 2	21.28 ± 7.41	19.92 ± 3.84	6.39	3.95	6.70	7.74
Media/Adven.	17.31 ± 6.72	16.18 ± 3.71	6.52	3.82	7.35	7.86

Discussion

- Vessel curvature and catheter eccentricity influence dose distribution of beta emitters
- Simplified models underestimate average doses as well as dose variability as compared to geometrically correct 3-D models
- No gold standard available
- However, tendencies shown in this study should prevail regardless of absolute values

Discussion

- Future extensions:
 - Plaque characterization
 - Monte-Carlo simulation
 - 4-D model by ECG sorting
 - ...

Conclusions

- Dose-delivery models for intravascular brachytherapy should consider vessel curvature and eccentricity of catheter
- Centering of the catheter of major interest to ensure prescribed dose is actually delivered homogeneously
- Our method may help improving both dosing models and delivery systems